# VanishPoint® Syringe

Automated Retraction Technology

## VanishPoint® Injektionsspritze

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



- The VanishPoint® syringe is a sterile, single-use syringe with an integrated retractable needle designed to reduce the risk of contaminated needlestick injuries.
- The intended use of the VanishPoint® syringe is to provide a safe and reliable method of preparing and administering injectable medication using a syringe and needle.
- Intended users of VanishPoint® syringes include healthcare professionals or laypersons, in healthcare settings or home settings.
- Indicated for administration of injectable medication dosages measured in milliliters (mL). Syringe graduation scale is not intended for insulin dose measurement or administration.
- Not indicated for administration of injectable insulin.
- Clinical benefits include accurate dose delivery of injectable medications and reduced risk of contaminated needlestick injuries. Activation of the needle retraction also reduces the risk of reuse and allows for safe, efficient disposal.
- May be used for any patient population with consideration given to patient size, prescribed medication dose, and route of administration.

### Instructions for Use:

1. Prepare and give injection using aseptic technique according to institutional policy.
2. Continue depressing plunger to activate automatic needle retraction **while needle is still in patient. Full dose is administered when the plunger handle is fully depressed**, which activates the needle retraction mechanism.
3. Needle will automatically retract into syringe, preventing exposure to contaminated needle and rendering syringe non-reusable. In the event that needle retraction mechanism does not activate, discard syringe in an appropriate sharps container per protocol of institution. Do not recap contaminated needles.
4. Dispose of VanishPoint® syringe in an appropriate sharps container per protocol of institution.

### Precautions:

- Single use only. Reuse of this device may result in exposure to bloodborne pathogens, including Hepatitis B virus (HBV), Hepatitis C virus (HCV), and human immunodeficiency virus (HIV).
- Contents are sterile, non-toxic, and non-pyrogenic. Do not use if product or package is damaged.
- Not made with natural rubber latex.
- Use only with attached needle. Needle cannot be changed.
- Automated needle retraction occurs only when barrel is emptied and plunger is fully depressed.
- For injections where full dose is not administered, expel remaining contents according to institutional policy if applicable and activate needle retraction.
- This device may be restricted for sale by or on the order of a physician.

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© 2021, 2024. Retractable Technologies, Inc. All rights reserved.