

Description	Needle Size	Units per Case	Boxes per Case	Catalog Number
VanishPoint IV Catheter - Radiopaque PUR	24G x 3/4" (19mm)	200	4	31221
VanishPoint IV Catheter - Radiopaque FEP*	24G x 3/4" (19mm)	200	4	35221
VanishPoint IV Catheter - Radiopaque PUR	22G x 1" (25mm)	200	4	31331
VanishPoint IV Catheter - Radiopaque FEP*	22G x 1" (25mm)	200	4	35331
VanishPoint IV Catheter - Radiopaque PUR	20G x 1 1/4" (32mm)	200	4	31441
VanishPoint IV Catheter - Radiopaque FEP*	20G x 1 1/4" (32mm)	200	4	35441
VanishPoint IV Catheter - Radiopaque PUR	18G x 1 1/4" (32mm)	200	4	31541
VanishPoint IV Catheter - Radiopaque FEP*	18G x 1 1/4" (32mm)	200	4	35541
VanishPoint IV Catheter - Radiopaque PUR	16G x 1 1/4" (32mm)	200	4	31641
VanishPoint IV Catheter - Radiopaque FEP*	16G x 1 1/4" (32mm)	200	4	35641
VanishPoint IV Catheter - Radiopaque PUR	14G x 1 1/4" (32mm)	200	4	31741
VanishPoint IV Catheter - Radiopaque FEP*	14G x 1 1/4" (32mm)	200	4	35741

*Contact your sales representative for current lead time

