



ADVANTAGES

- Accurate dose delivery
- Designed for patient comfort
- Easy to use
- Safe and effective
- Both hands remain behind the sharp
- Reduces risk of needlestick injury
- Non-reusable when activated
- Remains safe through disposal
- Minimizes disposal volume

**Saves Money,
Time, and Lives**



U-100 INSULIN SAFETY SYRINGES

| Description | Needle Size | # of Units per Case | # of Boxes per Case | Catalog Number |
|---|---------------------|---------------------|---------------------|----------------|
| 0.5mL VanishPoint U-100 Insulin Syringe (half-unit) | 30G x 3/16" (5mm) | 800 | 8 | 15391 |
| 0.5mL VanishPoint U-100 Insulin Syringe (half-unit) | 30G x 5/16" (8mm) | 800 | 8 | 15371 |
| 0.5mL VanishPoint U-100 Insulin Syringe | 30G x 3/16" (5mm) | 800 | 8 | 15291 |
| 0.5mL VanishPoint U-100 Insulin Syringe | 30G x 5/16" (8mm) | 800 | 8 | 15271 |
| 0.5mL VanishPoint U-100 Insulin Syringe | 30G x 1/2" (12.7mm) | 800 | 8 | 15221 |
| 1mL VanishPoint U-100 Insulin Syringe | 30G x 3/16" (5mm) | 800 | 8 | 10291 |
| 1mL VanishPoint U-100 Insulin Syringe | 30G x 5/16" (8mm) | 800 | 8 | 10271 |
| 1mL VanishPoint U-100 Insulin Syringe | 29G x 5/16" (8mm) | 800 | 8 | 10281 |
| 1mL VanishPoint U-100 Insulin Syringe | 29G x 1/2" (12.7mm) | 800 | 8 | 10211 |

Other available safety products



EasyPoint[®]
Needle

VanishPoint[®]
Blood Collection
Set

EasyPoint[®]
Blood Collection
Tube Holder
with Needle

VanishPoint[®]
IV Catheter

To order these products,
please contact
rtisales@retractable.com



RETRACTABLE
TECHNOLOGIES, INC.[®]

511 Lobo Lane • Little Elm • Texas • 75068-5295 • USA
Tel.: (972) 294-1770 • Fax: (972) 294-4400 • Toll Free: 1-888-703-1010

©2025 Retractable Technologies, Inc. All rights reserved. MKT-BRO-VPINS-0006 REV.0 EN, AUG 2025




**Safe and
accurate**



www.retractable.com



RETRACTABLE
TECHNOLOGIES, INC.[®]



Images based on actual sizes

Shorter needles may NOT require a pinch-up technique, reducing the risk of needlestick injury

NEEDLESTICK PREVENTION

- Safe and effective
- Reduces risk of needlestick injury
- Both hands remain behind the sharp

Real Safety[™]

Integrated safety mechanism allows for optimum view of syringe scale, promoting dose preparation ease & accuracy

VanishPoint[®] Insulin syringes are available with half unit increments for precise small dose delivery

Automated pre-removal retraction effectively reduces the risk of needlestick injuries.

VanishPoint[®] Insulin syringes are available with 5mm, 8mm and 12.7mm needles to ensure subcutaneous insulin delivery and maximize patient comfort

