

SAVE MONEY with the **NEW EasyPoint® Blood Collection Tube Holder with Needle**



FEATURES

- Pre-attached needle - no assembly required
- Automated retraction provides instant safety
- Allows for in-vein activation, which effectively reduces exposure to the contaminated needle
- Entire needle is completely covered when retracted, preventing exposure to the back-end needle
- Prevents needle removal and tube holder reuse

The point
that saves
you money



BENEFITS

- Single-use holder complies with safety regulations and is ready to use right out of package - **Saves Money**
- Improved efficiency by reducing preparation time - **Saves Time**
- Effectively reduces the risk of needlestick injuries - **Saves Lives**

EasyPoint® Blood Collection Tube Holder with Needle

Catalog No.	Needle Gauge	Needle Length
26011	22G	1¼"
26021	21G	1¼"



www.retractable.com

To order these products,
please contact
rtisales@retractable.com



RETRACTABLE
TECHNOLOGIES, INC.®