

OTHER SAFETY PRODUCTS FROM
RETRACTABLE TECHNOLOGIES

easypoint
**Blood Collection Tube
Holders with Needle**

The EasyPoint Blood Collection Tube Holder with Needle allows for in-vein activation, which effectively reduces exposure to the contaminated needle. The entire needle is covered when retracted.



**SAVES MONEY,
TIME & LIVES**

BLOOD COLLECTION TUBE HOLDERS WITH NEEDLE	
(100 per Case)	
26011	22G x 1 1/4"
26021	21G x 1 1/4"



RETRACTABLE
TECHNOLOGIES, INC.®

**THE
NEEDLE THAT
DISAPPEARS**



**P R O D U C T
G U I D E**

Through innovation, education, and the development of safe and reliable money-saving medical devices, Retractable Technologies, Inc. strives to be a catalyst in reducing the worldwide spread of infectious diseases.

Blood Collection Sets 50 Box per/200 per Case		
Catalog No.	Product Description	Needle Size
24011	7" (178mm) with Luer adapter	25G x 3/4" (19mm)
24021	7" (178mm) with Luer adapter	23G x 3/4" (19mm)
24031	7" (178mm) with Luer adapter	21G x 3/4" (19mm)
24111	12" (305mm) with Luer adapter	25G x 3/4" (19mm)
24121	12" (305mm) with Luer adapter	23G x 3/4" (19mm)
24131	12" (305mm) with Luer adapter	21G x 3/4" (19mm)
24311	12" (305mm) without Luer adapter	25G x 3/4" (19mm)
24321	12" (305mm) without Luer adapter	23G x 3/4" (19mm)
24331	12" (305mm) without Luer adapter	21G x 3/4" (19mm)



VanishPoint®
Blood Collection Sets

Reduces the risk of needlestick injuries during blood collection
Provides automated in-vein retraction and one-handed activation, effectively reducing the risk of needlestick injuries. Upon needle retraction, the tubing is clamped, reducing the risk of exposure to blood.




VanishPoint® IV Catheters
A safety IV catheter with our patented automated retraction

It is easy to use and allows for one-handed activation. Does not require additional components such as sliding sheaths, metal clips, or activation buttons.

Products are available in a variety of sizes.
For complete product lists,
visit our website at www.retractable.com

To order these products,
please contact
rtisales@retractable.com



**THE
GLOBAL
POINT**



www.retractable.com

511 Lobo Lane • Little Elm • Texas 75068-5295 • USA
Tel.: (972) 294-1770 • Fax: (972) 294-4400 • Toll Free: 1-888-703-1010



RETRACTABLE
TECHNOLOGIES, INC.®



VanishPoint®

easypoint



VanishPoint[®]
SYRINGES

RTI's VanishPoint Syringe is easy to use, requires no additional steps, and allows for single-handed, passive activation. After retraction, it requires less disposal space than other safety syringes and prevents disposal-related needlestick injuries. It is available in various sizes, and with a variety of needle gauges and lengths.

**SAVES
MONEY,
TIME &
LIVES**

easypoint[™]
RETRACTABLE NEEDLE

The EasyPoint Retractable Needle features automated pre-removal retraction that has made RTI products the industry standard. This product can be used with luer lock, luer slip, or pre-filled syringes.

AVAILABLE VANISHPOINT SYRINGES

Reference catalog numbers next to descriptions for easy ordering.

Catalog No.	Product Description	
VanishPoint® U-100 insulin syringes 100 per box/800 per case		
15391	0.5mL VanishPoint insulin syringe (half-unit)	30G x 3/16"
15371	0.5mL VanishPoint insulin syringe (half-unit)	30G x 5/16"
15291	0.5mL VanishPoint insulin syringe	30G x 3/16"
15271	0.5mL VanishPoint insulin syringe	30G x 5/16"
15221	0.5mL VanishPoint insulin syringe	30G x 1/2"
10291	1mL VanishPoint insulin syringe	30G x 3/16"
10271	1mL VanishPoint insulin syringe	30G x 5/16"
10281	1mL VanishPoint insulin syringe	29G x 5/16"
10211	1mL VanishPoint insulin syringe	29G x 1/2"
1mL VanishPoint® syringes 100 per box/800 per case		
10141	1mL VanishPoint syringe NEW	30G x 5/16"
10131	1mL VanishPoint syringe	27G x 1/2"
10151	1mL VanishPoint syringe	25G x 5/8"
10161	1mL VanishPoint syringe	25G x 1"
VanishPoint® Allergy Syringe Tray 25 per tray/1000 per case		
10134	1mL VanishPoint allergy syringes	27G x 1/2"
10154	1mL VanishPoint allergy syringes	25G x 5/8"

Catalog No.	Product Description	
3mL VanishPoint® syringes 100 per box/600 per case		
13021	3mL VanishPoint syringe	27G x 1 1/2"
10301	3mL VanishPoint syringe	25G x 5/8"
10391	3mL VanishPoint syringe	25G x 1"
13011	3mL VanishPoint syringe	25G x 1 1/2"
10311	3mL VanishPoint syringe	23G x 1"
10321	3mL VanishPoint syringe	23G x 1 1/2"
10331	3mL VanishPoint syringe	22G x 1"
10341	3mL VanishPoint syringe	22G x 1 1/2"
10351	3mL VanishPoint syringe	21G x 1"
10361	3mL VanishPoint syringe	21G x 1 1/2"
10371	3mL VanishPoint syringe	20G x 1"
10381	3mL VanishPoint syringe	20G x 1 1/2"
5mL VanishPoint® syringes 100 per box/600 per case		
10mL VanishPoint® syringes 100 per box/600 per case		


**THE
NEEDLE THAT
DISAPPEARS**




For complete catalog lists,
please scan the QR code or
visit our website at
www.retractable.com

AVAILABLE EASYPPOINT NEEDLES

Reference catalog numbers next to descriptions for easy ordering.

NEEDLES ONLY			
Catalog No.	Needle Size	Catalog No.	Needle Size
For syringes 3mL and smaller		For syringes 10mL and smaller	
82211	30G x 1/2"	85291	25G x 5/8"
82091	25G x 5/8"	85211	25G x 1"
82011	25G x 1"	85011	25G x 1 1/2"
82031	23G x 1"	85231	23G x 1"
		85021	22G x 1"
		85031	22G x 1 1/2"
		85041	21G x 1"
		85051	21G x 1 1/2"
		85061	20G x 1"
		85071	20G x 1 1/2"
		85081	18G x 1"
		85091	18G x 1 1/2"

WITH 3mL SYRINGES		
	Catalog No.	Needle Size
	82391	25G x 5/8"
	82311	25G x 1"
	85311	25G x 1 1/2"
	82331	23G x 1"
	85321	22G x 1"
	85331	22G x 1 1/2"
	85341	21G x 1"
	85351	21G x 1 1/2"
	85361	20G x 1"
	85371	20G x 1 1/2"
	85381	18G x 1"
	85391	18G x 1 1/2"




MONEY

Most syringes contain residual medication after use, and that medication is wasted. VanishPoint low dead space (LDS) syringes minimize medication wastage and may allow you to obtain extra doses from multi-dose vials, while providing precise and accurate dose delivery to patients. Activated VanishPoint syringes and EasyPoint needles require less disposal space than most other safety needles/syringes reducing waste management costs.



TIME

Automated retraction devices save time by streamlining workflow. By improving efficiency and safety, automated retraction saves valuable time.



LIVES

Retractable Technologies is dedicated to designing and producing products that create a safer environment for clinicians and patients. By using our safety products, healthcare facilities can substantially reduce the risk of needlestick injuries to employees, along with reducing costs.

