



1mL VanishPoint[®] Low Dead Space Safety Syringes

IMMUNIZATION **Superstar**

FEATURES

- Unique low dead space, attached needle design
- Integrated safety mechanism
- Automated needle retraction
- Pre-removal passive activation
- Safe and efficient disposal

BENEFITS

- As a result of medication savings, these syringes can result in an extra dose from a multi-dose vial.*
- Clinicians reported the ability to access extra doses from COVID-19 vaccine vials when using VanishPoint[®] 1mL Low Dead Space safety syringes.
- The integrated safety mechanism allows for a clear, unobstructed view of graduations, making dose preparation/administration easy and accurate.
- VanishPoint[®] syringes are automatically activated when the plunger handle is fully depressed – no additional steps or changes in hand position are required.
- Pre-removal activation effectively reduces exposure to the contaminated needle and the risk of needlestick injuries.
- VanishPoint[®] syringes have demonstrated superior safety and efficiency in acute care, long-term care, home care, and immunization programs.
- Activated VanishPoint[®] syringes are non-reusable, prevent disposal-related injuries, and significantly reduce disposal volume.

* Based on a multi-dose vial labeled for ten 0.5mL doses (0.05mL saved x 10 doses = 0.5mL, or one extra dose).



To order these products,
please contact
rtisales@retractable.com

