



- ▶ Automated needle retraction
- ▶ Effectively reduces the risk of needlestick injuries
- ▶ Quick and easy to use
- ▶ Activation can be achieved in a variety of ways
- ▶ May be activated before or after removal from patient
- ▶ Pre-removal activation reduces the risk of exposure to the contaminated sharp
- ▶ Activation of safety mechanism automatically covers syringe luer tip, reducing risk of exposure to medication, blood/body fluids
- ▶ Remains safe through disposal



compatible with both luer lock & luer slip syringes

may be used with pre-filled syringes



(50 per box /400 per case)

NEEDLES ONLY				WITH 3mL SYRINGES		
Catalog #	Needle Size	Catalog #	Needle Size		Catalog #	Needle Size
For syringes 3mL and smaller		For syringes 10mL and smaller			82391	25G x 5/8"
82211	30G x 1/2"	85291	25G x 5/8"		82311	25G x 1"
82091	25G x 5/8"	85211	25G x 1"		85311	25G x 1 1/2"
82011	25G x 1"	85011	25G x 1 1/2"		82331	23G x 1"
82031	23G x 1"	85231	23G x 1"		85321	22G x 1"
		85021	22G x 1"		85331	22G x 1 1/2"
		85031	22G x 1 1/2"		85341	21G x 1"
		85041	21G x 1"		85351	21G x 1 1/2"
		85051	21G x 1 1/2"		85361	20G x 1"
		85061	20G x 1"		85371	20G x 1 1/2"
		85071	20G x 1 1/2"		85381	18G x 1"
		85081	18G x 1"		85391	18G x 1 1/2"
		85091	18G x 1 1/2"			

OTHER AVAILABLE SAFETY PRODUCTS

- ▶ VanishPoint® Syringe
- ▶ VanishPoint® Blood Collection Tube Holder
- ▶ VanishPoint® Blood Collection Set
- ▶ VanishPoint® IV Catheter
- ▶ EasyPoint® Blood Collection Tube Holder with Needle



RETRACTABLE
TECHNOLOGIES, INC.®

To order these products, please contact
rtisales@retractable.com

511 Lobo Lane • Little Elm
Texas • 75068-5295 • USA
Tel.: (972) 294-1770
Fax: (972) 294-4400
Toll Free: 1-888-703-1010

© 2015, 2025 Retractable Technologies, Inc. All rights reserved. MKT-BR0-EPN-0002 REV.1 EN, NOV 2025



SAVES MONEY, TIME, AND LIVES



www.retractable.com




RETRACTABLE
TECHNOLOGIES, INC.®

easypoint™

RETRACTABLE NEEDLE

SAVES MONEY, TIME, AND LIVES



Available in a wide variety of
needle lengths & gauges



easy to use

- Minimal change in injection technique
- Allows for single-handed activation
- Clinician's hands and fingers remain behind the exposed sharp at all times

1



Securely attach EasyPoint needle to a syringe. Administer medication according to institutional policy.

2



Apply pressure to the color-coded tab and the safety chamber using a squeezing motion, allowing the needle to automatically retract into the safety chamber.

3



Dispose in an appropriate sharps container.