

VanishPoint® IV Catheters 50 per box/200 per case

Catalog Number	Gauge	Length	Catheter Material	Color Code
31221	24	¾"	Radiopaque PUR	Yellow
35221	24	¾"	Radiopaque FEP*	Yellow
31331	22	1"	Radiopaque PUR	Blue
35331	22	1"	Radiopaque FEP*	Blue
31441	20	1¼"	Radiopaque PUR	Pink
35441	20	1¼"	Radiopaque FEP*	Pink
31541	18	1¼"	Radiopaque PUR	Green
35541	18	1¼"	Radiopaque FEP*	Green
31641	16	1¼"	Radiopaque PUR	Grey
35641	16	1¼"	Radiopaque FEP*	Grey
31741	14	1¼"	Radiopaque PUR	Orange
35741	14	1¼"	Radiopaque FEP*	Orange

18-24 gauge IV catheters are suitable for use with power injectors up to 325 psi

*Special order items - please contact RTI for ordering information.

Other Available Safety Products



VanishPoint®
Syringe



EasyPoint®
Retractable
Needle



VanishPoint®
Blood Collection
Set



EasyPoint® Blood
Collection Tube
Holder with needle

- Integrated safety mechanism
- Minimizes exposure to the contaminated needle
- Allows for one-handed activation
- Once activated, the needle is safely retracted through disposal
- Easy to use

To order these products,
please contact
rtisales@retractable.com



RETRACTABLE
TECHNOLOGIES, INC.®

© 2012, 2025
Retractable Technologies, Inc.
All rights reserved MKT-BRO-
VPICV-0003 Rev.0 EN, AUG.2025

511 Lobo Lane, Little Elm
Texas 75068-5295 USA
Tel.: (972) 294-1010
Fax: (972) 294-4400
Toll Free: 1-888-703-1010

The New Standard for Safety



www.retractable.com

**SAVES MONEY,
TIME AND
LIVES**



RETRACTABLE
TECHNOLOGIES, INC.®

A safety
IV Catheter
featuring the
patented
VanishPoint[®]
automated
retraction!

SAFE AND EFFECTIVE

Sharp introducer needle
for insertion ease and
patient comfort

Push-off tab facilitates one handed
advancement and indicates bevel orientation

Clear finger grips and flashback
chamber allow for handling ease
and flashback visualization

Color coded on both
catheter hub and housing

Automatic needle retraction occurs with slight depression of the needle housing.
Once activated, the needle is safely retracted into the housing, reducing the risk of needlestick injuries

Perform venipuncture
observing for blood return.
Advance the catheter slightly
to assure full catheter entry
into the vein.

Holding the flashback
chamber stationary,
disengage the hub from
the introducer needle and
advance the catheter into
the vein.

Before withdrawing the
introducer needle from the
catheter, depress the color
coded end of the needle
housing to activate the
automated retraction.

The introducer needle will
retract into the housing.
Dispose of the introducer
needle in an appropriate
sharps container.